



**Xpod<sup>®</sup>**  
**Model 3011LP & 3012LP**  
**Specification and Technical Information**

**NONIN<sup>®</sup> Medical, Inc.**  
**13700 1st Avenue North**  
**Plymouth, Minnesota 55441-5443**  
**USA**

763-553-9968  
800-356-8874 (USA and Canada)  
Fax 763-553-7807  
E-mail: [info@nonin.com](mailto:info@nonin.com)  
[www.nonin.com](http://www.nonin.com)

© 2011 NONIN<sup>®</sup> Medical, Inc.

7602-000-03

# Table of Contents

<b>Specifications</b> .....	<b>1</b>
<b>Inputs/Outputs: Host Cable Conductors</b> .....	<b>4</b>
Inputs.....	4
Outputs.....	4
Data Format Selection .....	4
Select Data Format by Hardware Resistor .....	4
Select Data Format by Software Select.....	5
Patient Algorithm - SmartPoint™ .....	5
Serial Data Format #1 .....	6
Serial Data Format #2 .....	7
Serial Data Format #7 .....	11
Serial Data Format #8 .....	15
<b>Indications for Use</b> .....	<b>17</b>
<b>Contraindications</b> .....	<b>17</b>
<b>Warnings</b> .....	<b>17</b>
<b>Cautions</b> .....	<b>17</b>
<b>Accessories</b> .....	<b>19</b>
<b>Equipment Response Time</b> .....	<b>20</b>
<b>Testing Summary</b> .....	<b>21</b>
SpO <sub>2</sub> Accuracy Testing.....	21
Pulse Rate Motion Testing.....	21
Low Perfusion Testing .....	21
<b>Manufacturer's Declaration</b> .....	<b>22</b>



# Specifications

<b>Displayed Oxygen Saturation Range (SpO<sub>2</sub>)</b>	0 to 100%	
<b>Displayed Pulse Rate Range</b>	18 to 321 beats per minute (BPM)	
<b>Measurement Wavelengths and Output Power*</b>	Red: 660 nanometers @ 0.8 mW max. avg. Infrared: 910 nanometers @ 1.2 mW max. avg. (using Nonin PureLight <sup>®</sup> sensor)	
<b>SpO<sub>2</sub> Accuracy (A<sub>rms</sub>**)</b>	70 to 100%	
<b>No Motion</b>	<b>Adults/Pediatrics***</b>	<b>Neonates</b>
REUSABLE:		
8000AX Series:	± 2 digits	N/A
800XJ Series:	± 3 digits	N/A
8000SX Series:	± 2 digits	N/A
8000R:	± 3 digits	N/A
8000Q2:	± 3 digits	N/A
DISPOSABLE:		
6000 Series:	± 2 digits	± 3 digits
7000 Series:	± 2 digits	± 3 digits
6500 Series:	± 2 digits	N/A
<b>Motion</b>		
REUSABLE:		
8000AX Series:	± 2 digits	N/A
800XJ Series:	± 3 digits	N/A
8000SX Series:	± 3 digits	N/A
<b>Low Perfusion****</b>	± 2 digits	± 3 digits
<b>Pulse Rate Accuracy</b>	<b>Adults/Pediatrics***</b>	<b>Neonates</b>
<b>No Motion (18 – 300 BPM)</b>		
REUSABLE:		
8000AX Series:	± 3 digits	N/A
800XJ Series:	± 3 digits	N/A
8000SX Series:	± 3 digits	N/A
8000R:	± 3 digits	N/A
8000Q2:	± 3 digits	N/A
DISPOSABLE:		
6000 Series:	± 3 digits	± 3 digits
7000 Series:	± 3 digits	± 3 digits
6500 Series:	± 2 digits	N/A
<b>Motion (40 – 240 BPM)</b>		
REUSABLE:		
8000AX Series:	± 5 digits	N/A
800XJ Series:	± 5 digits	N/A
8000SX Series:	± 5 digits	N/A
<b>Low Perfusion (40 – 240 BPM) ****</b>	± 3 digits	± 3 digits

\* This information is especially useful for clinicians performing photodynamic therapy.

\*\*  $\pm 1$  Arms represents approximately 68% of measurements.

\*\*\* Includes infant patients.

\*\*\*\* Does not apply to those sensors listed as N/A under the neonate column, 8000R and 8000Q2.

**Notes:**

Reusable Group:

8000AX Finger Clip Sensors: 8000AA-1, 8000AA-3, 8000AP-1, 8000AP-3

8000XJ Flex Sensors: 8000J-1, 8000J-3, 8008J, 8001J

8000SX Soft Sensors: 8000SS, 8000SM, 8000SL

Disposable Group:

Flexi-Form<sup>®</sup> III (7000 Series) Sensors: 7000A, 7000P, 7000I, 7000N

6000C Series Sensors: 6000CA, 6000CP, 6000CI, 6000CN

Durafoam Disposable Sensors: 6500MA, 6500SA

**Temperature**

Operating: -05 to +50 °C for Xpod, 0 to +40 °C for sensor

Storage/Transportation: -40 to +70 °C

**Humidity**

Operating: 10% to 90% non-condensing

Storage/Transportation: 10% to 95% non-condensing

**Patient Isolation**

Type BF

**Leakage Current**

Not Applicable

**Dimensions**

Cable Length: 1 meter  $\pm 10\%$  or per customer request (2 meter maximum)

Housing: No larger than 1.1 cubic inches in volume +10%

**Weight**

No more than 75 grams

**Power Draw (typical)**

Typical 35mW or less with 3.3V input

Power Draw by Voltage Input:

Input Voltage	Power mW
1.0	50
1.2	46
1.4	41
1.6	39
1.8	37
2.0	36
2.2	36
2.4	35
2.6	34
2.8	34
3.0	33
3.2	33
3.3	34

Input Voltage	Power mW
3.4	34
3.6	35
3.8	36
4.0	37
4.2	39
4.4	40
4.6	41
4.8	43
5.0	45
5.2	46
5.4	47
5.5	47

**Inrush Power**

Shall not exceed an average of 275 mW within a sliding 5 mS window for any voltage input from 1.0 to 5.5 VDC.


**Voltage Input**

1.0 to 5.5 VDC, w/100 mV max. ripple

**Data I/O Signals**

0 to +3.7 VDC +/-5%

When connecting serial data directly to an MCU operating from 3.3V or when 3.3V is required, simply down convert 3.7V to 3.3V per suggested diagram:



**Patient Isolation**

Meets IEC 60601-1 Dielectric withstand  
Degree of protection: Type BF-Applied Part

**Leakage Current**

Not applicable

**Dimensions (approximate)**

53 mm (2.1 in.) x 20 mm (0.8 in.) x 15 mm (0.6 in.)

**Weight**

No greater than 75 g (2.7 oz.) including cable and Hirose connector

**Fluid Spill Resistance**

IP33

**Ruggedness**

Shock  
Ruggedness (Vibration)

IEC 60068-2-27  
Sinusoidal – IEC 60068-2-6  
Random – IEC 60068-2-64

**Sensors**

Designed to use Nonin-branded PureLight<sup>®</sup> sensors only (see Accessories)

**Shielding**

An RF shield is included (placed over the analog components)

**Regulatory Status/Information**

Contact Nonin [regulatory@nonin.com](mailto:regulatory@nonin.com)

**ROHS**

2002/95/EC

# Inputs/Outputs: Host Cable Conductors

## Inputs

- Red Wire = Power Input (1.0 to 5.5 VDC)
- Black Wire = Circuit Ground/Cable Shield
- Yellow Wire = Option 1 – Select data format using resistor in host device per the Data Format Selection.  
Option 2 – Select data format using software command per the Data Format Selection.

## Outputs

- Green Wire = Serial Output:  
9600 Baud, 8 data bits, One Start bit (Start bit =0), One Stop bit (Stop bit = 1), No Parity.  
Output Level: TTL (0 to 3.8 VDC)  
Host Input Impedence must greater or equal to 4kΩ.
- Blue Wire = PPG Output: Digital Pulse Indicator  
See Technical Note T-0604 for more information.
- Orange Wire = NC (reserved for future use)

---

**Note:** Xpod is not isolated from input voltage.

---

## Data Format Selection

The default data format for model 3011 is data format 1.

If you ordered the model 3011 and can not use the default DF1, the data format can be changed by using the hardware select resistor or software select data format options.

The default data format for model 3012 is data format 2.

If you ordered the model 3012 and can not use the default DF2, the data format can be changed by using the hardware select resistor or software select data format options.

For data formats 7 and 8, the host product must select the data format using the hardware or software data format select.

### Select Data Format by Hardware Resistor

Placing a resistor between Serial Input (yellow wire) and ground at the host connector end of cable can be used to select the data format. The table below lists the resistor values used per Data Format.

Output Format	Data Rate	Resistor Value (Ohms)
Data Format 1*	3 Bytes/once per second	Less than or equal to 626
Data Format 2	5 Bytes/75 times per second	8.2K +/- 5%
Data Format 7	5 Bytes/75 times per second	4.3K +/- 5%
Data Format 8	4 Bytes/once per second	22K +/- 5%

\*Data format 1 is retained for legacy purposes. If once per second data is desired, please use Data Format 8.

If the resistance is equal or greater to 297K and the software select option is not used, the default data format will be selected.

## Select Data Format by Software Select

If the host system does not use a resistor to select the data format, the host can select the Xpod data format by transmitting a 3-byte serial command within 1 second from power on of the Xpod. To select the data format, the host must send the Xpod the command as described below:

Byte #1: \$53 (ASCII value for letter capital "S"...for "Soft Select")

Byte #2: The Data Format (Hex value 01 for DF1, 02 for DF2, 07 for DF7, 08 for DF8)

Byte #3: Checksum (Hex value) = Byte#1 + Byte#2

The Xpod must receive the user command within 1 second after power is applied to the Xpod. Commands after the first second from power on will not be processed.

Data format select by serial configuration from host to XPOD:

9600 Baud, 8 data bits, One Start bit (Start bit =0), One Stop bit (Stop bit = 1), No Parity.

Key Points when using the software select feature:

1. The data format selected by the software command will be lost when power to the Xpod is removed.
2. When using the software select feature, make sure the host device sends the command within one second after applying power to the Xpod.
3. If the host receives the wrong data format, remove power from the Xpod. Then apply power and send the desired software command within one second.

## Patient Algorithm - SmartPoint™

Data formats 2, 7, and 8 provide a SmartPoint indicator. The SmartPoint Algorithm qualifies the data for recording purposes and eliminates the guesswork of determining when the patient measurement is qualified for recording purposes. When the SmartPoint Algorithm indicates the reading is high quality, the SPA bit will be set in data formats 2, 7, and 8.

## Serial Data Format #1

This data format provides continuous data transmission of a 3 byte data packet sent once per second. The data packet includes real-time data including: SpO<sub>2</sub> and Pulse Rate formatted for display, and status of the measurement.

### Packet Description

Three bytes of data are transmitted 1 once per second.

Byte 1 - Status							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
1	SNSD	OOT	LPRF	MPRF	ARTF	HR8	HR7
<b>*Note:</b> Bit 7 is always set							

Byte 2 - Heart Rate							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	HR6	HR5	HR4	HR3	HR2	HR1	HR0
<b>*Note:</b> Bit 7 is always clear							

Byte 3 - SpO <sub>2</sub>							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	SP6	SP5	SP4	SP3	SP2	SP1	SP0
<b>*Note:</b> Bit 7 is always clear							

The following are all active high:

SNSD:	Sensor Disconnect	Sensor is not connected to oximeter or sensor is inoperable.
OOT:	Out Of Track	An absence of consecutive good pulse signals.
LPRF:	Low Perfusion	Amplitude representation of low/no signal quality.
MPRF:	Marginal Perfusion	Amplitude representation of low/marginal signal quality.
ARTF:	Artifact	Indicated artifact condition on each pulse.
HR8 – HR0:	Heart Rate	Standard 4-beat average values without display holds.
SP6 – SP0:	SpO <sub>2</sub>	Standard 4-beat average values without display holds.

These SpO<sub>2</sub> and HR values are formatted for recording purposes and are updated every 1/3 of second. When the sensor is removed from the site, these values will be formatted with the missing data value. The following output options are available in standard mode:

HR:                   4-beat Pulse Rate Average  
 SpO<sub>2</sub>:               4-beat SpO<sub>2</sub> Average

When SpO<sub>2</sub> and HR values cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127. The missing data could be result of these conditions:

1. Sensor is positioned improperly.
2. Sensor was removed prior to a reading.
3. Signal at the sensor site is not discernable. Warm the site or choose a different site.

## Serial Data Format #2

This data format provides continuous data transmission of a 5 byte data packet sent 75 times per second. The data packet includes real-time data including: 8-bit waveform value, six different output options for the SpO<sub>2</sub> value, four different averaging options for the Pulse Rate values, and options formatted for both recording and display purposes, as well as status information for the measurement.

### Packet Description

A frame consists of 5 bytes; a packet consists of 25 frames. Three packets (75 frames) are transmitted each second.

Packet	Frame				
	Byte 1	Byte 2	Byte 3	Byte 4	Byte 5
1	01	STATUS	PLETH	HR MSB	CHK
2	01	STATUS	PLETH	HR LSB	CHK
3	01	STATUS	PLETH	SpO <sub>2</sub>	CHK
4	01	STATUS	PLETH	SREV	CHK
5	01	STATUS	PLETH	reserved	CHK
6	01	STATUS	PLETH	reserved	CHK
7	01	STATUS	PLETH	reserved	CHK
8	01	STATUS	PLETH	STAT2	CHK
9	01	STATUS	PLETH	SpO <sub>2</sub> -D	CHK
10	01	STATUS	PLETH	SpO <sub>2</sub> Fast	CHK
11	01	STATUS	PLETH	SpO <sub>2</sub> B-B	CHK
12	01	STATUS	PLETH	reserved	CHK
13	01	STATUS	PLETH	reserved	CHK
14	01	STATUS	PLETH	E-HR MSB	CHK
15	01	STATUS	PLETH	E-HR LSB	CHK
16	01	STATUS	PLETH	E-SpO <sub>2</sub>	CHK
17	01	STATUS	PLETH	E-SpO <sub>2</sub> -D	CHK
18	01	STATUS	PLETH	reserved	CHK
19	01	STATUS	PLETH	reserved	CHK
20	01	STATUS	PLETH	HR-D MSB	CHK
21	01	STATUS	PLETH	HR-D LSB	CHK
22	01	STATUS	PLETH	E-HR-D MSB	CHK
23	01	STATUS	PLETH	E-HR-D LSB	CHK
24	01	STATUS	PLETH	reserved	CHK
25	01	STATUS	PLETH	reserved	CHK

#### Notes:

Byte number 1 in each frame is set to a value of 1.

Reserved bytes are undefined (range of 0 to 255).

**Byte 1 – START BYTE:**

Always set to a 01 value.

**Byte 2 – STATUS BYTE**

This byte provides status information at a rate of 1/75 of a second.

Range: 128 to 255

Byte 2 - Status							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
1	SNSD	ARTF	OOT	SNSF	YPRF		SYNC
					RPRF	GPRF	
*Note: Bit 7 is always set.							

The following are all active high:

SNSD:	Sensor Disconnect	Sensor is not connected to oximeter or sensor is inoperable.
ARTF:	Artifact – short term	Indicates artifact condition of each pulse (occurs only during pulse).
OOT:	Out Of Track	An absence of consecutive good pulse signals.
SNSA:	Sensor Alarm	Device is providing unusable data for analysis ( set when the finger is removed or sensor is disconnected).
RPRF:	*Red Perfusion	Amplitude representation of low/poor signal quality (occurs only during pulse).
YPRF:	*Yellow Perfusion	Amplitude representation of low/marginal signal quality (occurs only during pulse).
GPRF:	*Green Perfusion	Amplitude representation of high signal quality (occurs only during pulse).
SYNC:	Frame Sync	1 on Frame 1 (0 on frames 2 through 25).

\* The oximeter reports each pulse by setting/clearing the RPRF and GPRF bits for a period of 12 frames (160 ms). The table below describes the condition and state of the pulse perfusion bits.

Condition	RPRF Bit 2 of Status Byte	GPRF Bit 1 of Status Byte
Green – high pulse signal	0	1
Yellow – low/marginal pulse signal	1	1
Red – low/no pulse signal	1	0

### Byte 3 – PLETH BYTE

This byte consists of an 8 bit plethysmographic waveform (pulse waveform). The pulse oximeter infra-red signal is filtered and then compressed into an 8 bit value. The compression provides good detail for low to medium pulse signals. For an uncompressed waveform with better resolution, refer to Data Format 7.

Range: 0 to 255

### Byte 4 – FLAT BYTE

This byte is used for SpO<sub>2</sub>, Pulse Rate, and information that can be processed at a rate of 1/3 of a second.

Range: 00 to 127

SREV: Oximeter Firmware Revision Level

STAT2: Status Byte 2 (occurs 1 of 25) - description given below

Byte 4 – STAT 2							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	R	SPA	R	R	R	R	R
*Note: Bit 7 is always set.							

The following are all active high:

SPA: High Quality SmartPoint Measurement  
 R: Reserved (range - 0 or 1), for future use

### Standard Mode - Formatted for Recording Purposes:

These values are formatted for recording purposes and are updated every 1/3 of second. When the finger is removed from the device these values will be formatted with the missing data value. The following output options are available in standard mode:

HR: 4-beat Pulse Rate Average  
 E-HR: 8-beat Pulse Rate Extended Average  
 SpO<sub>2</sub>: 4-beat SpO<sub>2</sub> Average  
 E-SpO<sub>2</sub>: 8-beat SpO<sub>2</sub> Extended Average  
 SpO<sub>2</sub> Fast: 4-beat Average optimized for fast responding  
 SpO<sub>2</sub> B-B: Beat to Beat value – No Average

When SpO<sub>2</sub> and HR cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127.

*Display Mode - Formatted for Display Purposes:*

These values are formatted for display purposes and are updated every 1.5 seconds. When the sensor is removed from the site, the last SpO<sub>2</sub> and Pulse Rate reading will be reported for 10 seconds before changing to the missing data value. During this 10 second period the sensor alarm bit (SNSA) is set, indicating that the sensor has been removed. This feature is useful for spot-check measurements. The following output options are available in Display Mode:

- HR-D: 4-beat Pulse Rate Average
- E-HR-D: 8-beat Pulse Rate Extended Average
- SpO<sub>2</sub>-D: 4-beat SpO<sub>2</sub> Average
- E-SpO<sub>2</sub>-D: 8-beat SpO<sub>2</sub> Extended Average

When SpO<sub>2</sub> and HR cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127. The missing data could be result of these conditions:

1. Sensor is positioned improperly.
2. Sensor was removed prior to a reading.

Signal at the sensor site is not discernable. Warm the site or choose a different site.

HR Format:

	7	6	5	4	3	2	1	0
<b>HR MSB</b>	0	R	R	R	R	R	HR8	HR7

	7	6	5	4	3	2	1	0
<b>HR LSB</b>	0	HR6	HR5	HR4	HR3	HR2	HR1	HR0

SpO<sub>2</sub> Format:

	7	6	5	4	3	2	1	0
<b>SpO<sub>2</sub></b>	0	SP6	SP5	SP4	SP3	SP2	SP1	SP0

R = Reserved (range 0 or 1)

**Byte 5 – CHK**

This byte is used for the checksum of bytes 1 through 4.

Range: 00 to 255

$$\text{CHK: Checksum} = (\text{Byte 1}) + (\text{Byte 2}) + (\text{Byte 3}) + (\text{Byte 4}) \text{ modulo } 256$$

## Serial Data Format #7

This data format provides the same information as Data Format 2, except that the waveform value provides the full resolution of 16 bits instead of 8 bits. This data format must be selected by hardware resistor select or software select.

### Packet Description

A frame consists of 5 bytes; a packet consists of 25 frames. Three packets (75 frames) are transmitted each second.

Packet	Frame				
	Byte 1	Byte 2	Byte 3	Byte 4	Byte 5
1	STATUS	PLETH MSB	PLETH LSB	HR MSB	CHK
2	STATUS	PLETH MSB	PLETH LSB	HR LSB	CHK
3	STATUS	PLETH MSB	PLETH LSB	SpO <sub>2</sub>	CHK
4	STATUS	PLETH MSB	PLETH LSB	SREV	CHK
5	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
6	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
7	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
8	STATUS	PLETH MSB	PLETH LSB	STAT2	CHK
9	STATUS	PLETH MSB	PLETH LSB	SpO <sub>2</sub> -D	CHK
10	STATUS	PLETH MSB	PLETH LSB	SpO <sub>2</sub> Fast	CHK
11	STATUS	PLETH MSB	PLETH LSB	SpO <sub>2</sub> B-B	CHK
12	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
13	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
14	STATUS	PLETH MSB	PLETH LSB	E-HR MSB	CHK
15	STATUS	PLETH MSB	PLETH LSB	E-HR LSB	CHK
16	STATUS	PLETH MSB	PLETH LSB	E-SpO <sub>2</sub>	CHK
17	STATUS	PLETH MSB	PLETH LSB	E-SpO <sub>2</sub> -D	CHK
18	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
19	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
20	STATUS	PLETH MSB	PLETH LSB	HR-D MSB	CHK
21	STATUS	PLETH MSB	PLETH LSB	HR-D LSB	CHK
22	STATUS	PLETH MSB	PLETH LSB	E-HR-D MSB	CHK
23	STATUS	PLETH MSB	PLETH LSB	E-HR-D LSB	CHK
24	STATUS	PLETH MSB	PLETH LSB	reserved	CHK
25	STATUS	PLETH MSB	PLETH LSB	reserved	CHK

---

#### Notes:

Byte number 1 in each frame is greater than 127.

Reserved bytes are undefined (range of 0 to 255).

---

## Byte 1 – STATUS BYTE

This byte provides status information at a rate of 1/75 of a second.

Range: 128 to 255

Byte 1 - Status							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
1	SNSD	ARTF	OOT	SNSF	YPRF		SYNC
					RPRF	GPRF	
*Note: Bit 7 is always set.							

The following are all active high:

SNSD:	Sensor Disconnect	Sensor is not connected to oximeter or sensor is inoperable.
ARTF:	Artifact	Indicates artifact condition of each pulse (occurs only during pulse).
OOT:	Out Of Track	An absence of consecutive good pulse signals.
SNSA:	Sensor Alarm	Device is providing unusable data for analysis (set when the finger is removed or sensor is disconnected).
RPRF:	*Red Perfusion	Amplitude representation of low/no pulse signal (occurs only during pulse).
YPRF:	*Yellow Perfusion	Amplitude representation of low/marginal signal quality (occurs only during pulse).
GPRF:	*Green Perfusion	Amplitude representation of high signal quality (occurs only during pulse).
SYNC:	Frame Sync	= 1 to Frame 1 (=0 on frames 2 through 25).

\* The oximeter reports each pulse by setting/clearing the RPRF and GPRF bits for a period of 12 frames (160 ms). The table below describes the condition and state of the pulse perfusion bits.

Condition	RPRF Bit 2 of Status Byte	GPRF Bit 1 of Status Byte
Green – high pulse signal	0	1
Yellow – low/marginal pulse signal	1	1
Red – low/no pulse signal	1	0

### Byte 2 & 3 – PLETH BYTE

These two bytes consist of a 16 bit plethysmographic waveform (pulse waveform).

Range: 0 to 65535 (MSB:LSB )

Byte 2 = MSB Pulse Waveform

Byte 3 = LSB Pulse Waveform

Pulse waveform value = (Byte 2 decimal value \* 256) + Byte 3 decimal value

### Byte 4 – FLOAT BYTE

This byte is used for SpO<sub>2</sub>, Pulse Rate, and information that can be processed at a rate of 1/3 of a second.

Range: 00 to 127

SREV: Oximeter Firmware Revision Level  
 STAT2: Status Byte 2 (occurs 1 of 25) - description given below

Byte 4 – STAT 2							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	R	SPA	R	R	R	R	R

The following are all active high:

SPA: High quality SmartPoint Measurement  
 R: Reserved (range - 0 or 1), for future use

#### Standard Mode - Formatted for Recording Purposes:

These values are formatted for recording purposes and are updated every 1/3 of second. When the sensor is removed from the site, these values will be formatted with the missing data value. The following output options are available in standard mode:

HR: 4-beat Pulse Rate Average  
 E-HR: 8-beat Pulse Rate Extended Average  
 SpO<sub>2</sub>: 4-beat SpO<sub>2</sub> Average  
 E- SpO<sub>2</sub>: 8-beat SpO<sub>2</sub> Extended Average  
 SpO<sub>2</sub> Fast: 4-beat Average optimized for fast responding  
 SpO<sub>2</sub> B-B: Beat to Beat value – No Average

When SpO<sub>2</sub> and HR cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127. The missing data could be result of these conditions:

1. Sensor is positioned improperly.
2. Sensor was removed prior to a reading.
3. Signal at the sensor site is not discernable. Warm the site or choose a different site.

*Display Mode - Formatted for Display Purposes:*

These values are formatted for display purposes and are updated every 1.5 seconds. When the sensor is removed from the site, the last SpO<sub>2</sub> and Pulse Rate reading will be reported for 10 seconds before changing to the missing data value. During this 10 second period the sensor alarm bit (SNSA) is set, indicating that the sensor has been removed. This feature is useful for spot-check measurements. The following output options are available in Display Mode:

- HR-D: 4-beat Pulse Rate Average
- E-HR-D: 8-beat Pulse Rate Extended Average
- SpO<sub>2</sub>-D: 4-beat SpO<sub>2</sub> Average
- E- SpO<sub>2</sub>-D: 8-beat SpO<sub>2</sub> Extended Average

When SpO<sub>2</sub> and HR cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127.

HR Format:

	7	6	5	4	3	2	1	0
<b>HR MSB</b>	0	R	R	R	R	R	HR8	HR7

	7	6	5	4	3	2	1	0
<b>HR LSB</b>	0	HR6	HR5	HR4	HR3	HR2	HR1	HR0

SpO<sub>2</sub> Format:

	7	6	5	4	3	2	1	0
<b>SpO<sub>2</sub></b>	0	SP6	SP5	SP4	SP3	SP2	SP1	SP0

R: Reserved (range- 0 or 1)

**Byte 5 – CHK**

This byte is used for the checksum of bytes 1 through 4.

Range: 00 to 255

CHK: Checksum = (Byte 1) + (Byte 2) + (Byte 3) + (Byte 4) modulo 256

## Serial Data Format #8

This data format provides continuous data transmission of a 4 byte data packet sent once per second. The data packet includes real-time data including: SpO<sub>2</sub> and Pulse Rate formatted for display and status of the measurement. This data format must be selected by the hardware resistor select or software select feature.

### Packet Description

Three bytes of data are transmitted once per second.

Byte 1 - Status							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
1	SNSD	OOT	LPRF	MPRF	ARTF	HR8	HR7
<b>*Note:</b> Bit 7 is always set							

Byte 2 - Heart Rate (HR-D)							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	HR6	HR5	HR4	HR3	HR2	HR1	HR0
<b>*Note:</b> Bit 7 is always clear							

Byte 3 - SpO <sub>2</sub> -D							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	SP6	SP5	SP4	SP3	SP2	SP1	SP0
<b>*Note:</b> Bit 7 is always clear							

Byte 4 – Status2							
BIT7	BIT6	BIT5	BIT4	BIT3	BIT2	BIT1	BIT0
0	R	SPA	R	SNSA	R	R	R
<b>*Note:</b> Bit 7 is always clear							

The following are all active high:

SNSD:	Sensor Disconnect	Sensor is not connected to oximeter or sensor is inoperable.
ARTF:	Artifact	Indicated artifact condition on each pulse.
OOT:	Out Of Track	An absence of consecutive good pulse signals.
LPRF:	Low Perfusion	Amplitude representation of low/no signal quality.
MPRF:	Marginal Perfusion	Amplitude representation of low/marginal signal quality.
SNSA:	Sensor Alarm	Device is providing unusable data for analysis or sensor is disconnected.
SPA:	SmartPoint Algorithm	High quality SmartPoint measurement.
HR8 – HR0:	Heart Rate (HR-D)	4-beat Pulse Rate average formatted for display.
SP6 – SP0:	SpO <sub>2</sub> (SpO <sub>2</sub> -D)	4-beat SpO <sub>2</sub> average formatted for display.
R	Reserved (range – 0 or 1)	Reserved for future use.

The SpO<sub>2</sub> and Pulse Rate values are formatted for display purposes and are updated every 1.5 seconds. When the sensor is removed from the site, the last SpO<sub>2</sub> and Pulse Rate reading will be reported for 10 seconds before changing to the missing data value. During this 10 second period the sensor alarm bit (SNSA) is set, indicating that the sensor has been removed. This feature is useful for spot-check measurements. The following output options are available in Display Mode:

HR-D:	4-beat Pulse Rate Average
SpO <sub>2</sub> -D:	4-beat SpO <sub>2</sub> Average

When SpO<sub>2</sub> and HR cannot be computed, the system will send a missing data indicator. For missing data, the HR equals 511 and the SpO<sub>2</sub> equals 127. The missing data could be result of these conditions:

1. Sensor is positioned improperly.
2. Sensor was removed prior to a reading.
3. Signal at the sensor site is not discernable. Warm the site or choose a different site.

## Indications for Use

The Xpod is intended to provide medical device manufacturers with a small, low-power oximeter that can be easily attached to a host device externally. The Xpod measures functional oxygen saturation of arterial hemoglobin (%SpO<sub>2</sub>) and pulse rate (BPM) for adult, pediatric, infant and neonatal patients. When integrated with a medical device manufacturer's host system, the Xpod may be used in any environment where pulse oximetry measurements are made.

### Contraindications

- Do not use this device in an MR environment.
- Explosive Hazard: Do not use this device in an explosive atmosphere or in the presence of flammable anesthetics or gases
- This module does not meet defibrillation-proof requirement per IEC 60601-1: 1988/A1:1995, clause 17.h.

### Warnings

- Use only with Nonin-branded PureLight® pulse oximeter sensors. These sensors are manufactured to meet the accuracy specifications for Nonin pulse oximeters. Using other manufacturers' sensors can result in inaccurate pulse oximeter performance.
- Loss of monitoring can result if any objects hinder the pulse measurement. Ensure that no blood flow restrictors (e.g., blood pressure cuff) hinder pulse measurements.
- As with all medical equipment, carefully route cables and connections to reduce the possibility of entanglement or strangulation.
- Operation of this module below the minimum amplitude of 0.3% modulation may cause inaccurate results.
- The use of accessories, sensors, and cables other than those specified by Nonin may result in increased emission and/or decreased immunity of this device.
- Do not use a damaged sensor.

### Cautions

- The accuracy of the SpO<sub>2</sub> measurement may be affected if the total sensor cable length (including extension cables) is greater than 3 meters (9.8 feet).
- Follow local, state, or national governing ordinances and recycling instructions regarding disposal or recycling of the device and device components.
- In compliance with the European Directive on Waste Electrical and Electronic Equipment (WEEE) 2002/96/EC, do not dispose of this product as unsorted municipal waste. This device contains WEEE materials; please contact your distributor regarding take-back or recycling of the device.
- This pulse oximeter module is designed to determine the percentage of arterial oxygen saturation of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:
  - excessive ambient light
  - excessive motion
  - electro-surgical interference
  - blood flow restrictors (arterial catheters, blood pressure cuffs, infusion lines, etc.)
  - moisture in the sensor
  - improperly applied sensor
  - incorrect sensor type
  - poor pulse quality
  - venous pulsations
  - anemia or low hemoglobin concentrations
  - cardiogreen or other intravascular dyes
  - carboxyhemoglobin
  - methemoglobin

## Cautions

- dysfunctional hemoglobin
  - artificial nails or fingernail polish
  - a sensor not at heart level
- 
- This device has motion tolerant software that minimizes the likelihood of motion artifact being misinterpreted as good pulse quality. In some circumstances, however, this device may still interpret motion as good pulse quality. This covers all available outputs (i.e. SpO<sub>2</sub>, HR, PLETH, PPG).
- 
- Inspect the sensor application site at least every 6 to 8 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity may vary due to medical status or skin condition. Discontinue use of adhesive tape strips if the patient exhibits an allergic reaction to the adhesive material.
- 
- This equipment complies with IEC 60601-1-2:2004 for electromagnetic compatibility for medical electrical equipment and/or systems. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device. Medical electrical equipment needs special precautions regarding EMC, and all equipment must be installed and put into service according to the EMC information specified in this manual.
- 
- Portable and mobile RF communications equipment may affect medical electrical equipment.
- 
- Oximeter readings may be affected by the use of an electrosurgical unit (ESU)
- 
- The oximeter sensor may not work on cold extremities due to reduced circulation. Warm or rub the finger to increase circulation, or reposition the sensor.
- 
- A functional tester cannot be used to assess the accuracy of a pulse oximeter or sensor.

For more information about required safety and regulatory requirements for pulse oximeters, refer to ISO 9919: 2005 and IEC 60601-1: 1988/A1:1995. Additional safety information can be found in the labeling provided with each Nonin sensor.

## Accessories

The following Nonin accessories may be used with the Xpod module. See the respective sensor instructions for detailed information regarding specified sensor use (patient population, body/tissue, and application).

<b>Model Number</b>	<b>Description</b>
8000AA-1	Adult Articulated Internal Spring Finger Clip, 1 m (3 ft) cable
8000AA-3	Adult Articulated Internal Spring Finger Clip, 3 m (9.8 ft) cable
8000AP-1	Pediatric External Spring Finger Clip, 1 m (3 ft) cable
8000AP-3	Pediatric External Spring Finger Clip, 3 m (9.8 ft) cable
8000J-1	Adult Flex, 1 m (3 ft) cable
8000J-3	Adult Flex, 3 m (9.8 ft) cable
8001J	Neonatal Flex, 1 m (3 ft) cable
8008J	Infant Flex, 1 m (3 ft) cable
8000Q2	Ear Clip, 1 m (3 ft) cable
8000R	Reflectance, 1 m (3 ft) cable
8000SS	Sensor, Reusable, Soft, Small, 1 m (3 ft) cable
8000SM	Sensor, Reusable, Soft, Medium, 1 m (3 ft) cable
8000SL	Sensor, Reusable, Soft, Large, 1 m (3 ft) cable
7000A	Flexi-Form® II Adult, 1 m (3 ft) cable, 10-pack
7000P	Flexi-Form® II Pediatric, 1 m (3 ft) cable, 10-pack
7000I	Flexi-Form® II Infant, 1 m (3 ft) cable, 10-pack
7000N	Flexi-Form® II Adult, 1 m (3 ft) cable, 10-pack
6000CA	Sensor, Disposable, Adult, 45 cm (17.5 in) cable
6000CP	Sensor, Disposable, Pediatric, 45 cm (17.5 in) cable
6000CI	Sensor, Disposable, Infant, 90 cm (35.5 in) cable
6000CN	Sensor, Disposable, Neonate, 90 cm (35.5 in) cable
6500SA	Sensor, Durafoam Disposable, Standard, 1 m (3 ft) cable
6500MA	Sensor, Durafoam Disposable, Small, 1 m (3 ft) cable
UNI-RA-0	7.5" 90-degree Patient Cable
UNI EXT-X	Patient Extension Cable (select 1, 3, 6, or 9 meters)

## Equipment Response Time


SpO <sub>2</sub> Values	Average	Latency
Standard/Fast Averaged SpO <sub>2</sub>	4 beat average	2 beats
Extended Averaged SpO <sub>2</sub>	8 beat average	2 beats

Pulse Rate Values	Average	Latency
Standard/Fast Averaged Pulse Rate	4 beat average	2 beats
Extended Averaged Pulse Rate	8 beat average	2 beats

### Example – SpO<sub>2</sub> Exponential Averaging

SpO<sub>2</sub> decreases 0.75% per second (7.5% over 10 seconds)

Pulse Rate – 75 BPM



Specific to this example:

- The response of the 4-beat average is 1.5 seconds.
- The response of the 8-beat average is 3 seconds.

## Testing Summary

SpO<sub>2</sub> accuracy, motion and low perfusion testing was conducted by Nonin Medical, Incorporated as described below.

### SpO<sub>2</sub> Accuracy Testing

SpO<sub>2</sub> accuracy testing is conducted during induced hypoxia studies on healthy, non-smoking, light-to-dark-skinned subjects during motion and no-motion conditions in an independent research laboratory. The measured arterial hemoglobin saturation value (SpO<sub>2</sub>) of the sensors is compared to arterial hemoglobin oxygen (SaO<sub>2</sub>) value, determined from blood samples with a laboratory co-oximeter. The accuracy of the sensors in comparison to the co-oximeter samples measured over the SpO<sub>2</sub> range of 70 – 100%. Accuracy data is calculated using the root-mean-squared ( $A_{rms}$  value) for all subjects, per ISO 9919:2005, Standard Specification for Pulse Oximeters for Accuracy.

### Pulse Rate Motion Testing

This test measures pulse rate accuracy with motion artifact simulation introduced by a pulse oximeter tester. This test determines whether the oximeter meets the criteria of ISO 9919:2005 for pulse rate during simulated movement, tremor, and spike motions.

### Low Perfusion Testing

This test uses an SpO<sub>2</sub> Simulator to provide a simulated pulse rate, with adjustable amplitude settings at various SpO<sub>2</sub> levels. The module must maintain accuracy in accordance with ISO 9919:2005 for pulse rate and SpO<sub>2</sub> at the lowest obtainable pulse amplitude (0.3% modulation).

## Manufacturer's Declaration

See the following tables for specific information regarding this module's compliance to IEC 60601-1-2:2007.


**Table 1: Electromagnetic Emissions**

Emissions Test	Compliance	Electromagnetic Environment—Guidance
<i>This module is intended for use in the electromagnetic environment specified below. The customer and/or user of this device should ensure that it is used in such an environment.</i>		
RF Emissions CISPR 11	Group 1	This module uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	This module is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	N/A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	N/A	

**Table 2: Electromagnetic Immunity**

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment—Guidance
<i>This module is intended for use in the electromagnetic environment specified below. The customer and/or user of this device should ensure that it is used in such an environment.</i>			
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	N/A	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	N/A	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	N/A	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the module requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery pack.
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>Note:</b> $U_T$ is the AC mains voltage before application of the test level.			

**Table 3: Guidance and Manufacturer’s Declaration—Electromagnetic Immunity**

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment—Guidance
<p><i>This module is intended for use in the electromagnetic environment specified below. The customer and/or user of this module should ensure that it is used in such an environment.</i></p>			
<p>Portable and mobile RF communications equipment should be used no closer to any part of the module, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>			
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p><b>Recommended Separation Distance</b></p> <p><math>d = 1.17 \sqrt{P}</math></p> <p><math>d = 1.17 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.33 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (<math>W</math>) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p><b>Notes:</b></p> <ul style="list-style-type: none"> <li>• At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</li> <li>• These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</li> </ul> <p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the module.</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [3] V/m.</p>			

**Table 4: Recommended Separation Distances**

The following table describes the recommended separation distances between portable and mobile RF communications equipment and this module.

<p><i>This module is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Customers or users of this module can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the module as recommended below, according to maximum output power of the communications equipment.</i></p>			
<p>Rated Maximum Output Power of Transmitter</p> <p>W</p>	<p>Separation Distance According to Frequency of Transmitter</p>		
	<p>150 kHz to 80 MHz</p> <p><math>d = 1.17\sqrt{P}</math></p>	<p>80 MHz to 800 MHz</p> <p><math>d = 1.17\sqrt{P}</math></p>	<p>800 MHz to 2.5 GHz</p> <p><math>d = 2.33\sqrt{P}</math></p>
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance <math>d</math> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p>			
<p><b>Notes:</b></p> <ul style="list-style-type: none"> <li>• At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</li> <li>• These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</li> </ul>			