

## 2500 FAQ'S

### **Q. What is included in the 2500 package assembly?**

- A.**
- 2500 Pulse Oximeter • 2500 Operator's Manual on CD
  - 8000AA Pulse Oximeter Sensor • 4 AA-Size Alkaline Batteries

### **Q. How do you know you're getting an accurate reading?**

**A.** The tricolor LED display (m) provides a visual indication of the pulse signal quality, while blinking at the corresponding pulse rate. This display changes colors to alert you to changes in pulse quality that may affect the readings: *green* indicates a good pulse signal, *yellow* indicates a marginal pulse signal, and *red* indicates an inadequate pulse signal. Allow the 2500A to stabilize, observing about 10 seconds of continuous green-colored pulse quality before relying on the displayed values. It is common for the displayed values to fluctuate slightly over a period of several seconds. If the pulse quality display blinks yellow or red, try another finger

### **Q. What kind of battery is used in the 2500?**

**A.** Four AA Alkaline or you can purchase the 2500C-UNIV Rechargeable NiMH Battery System.

### **Q. How many hours of battery life can I expect?**

**A.** Four AA alkaline batteries provide the device with approximately 80 hours of continuous operation, while the Rechargeable NiMH Battery Pack provides approximately 40 hours of continuous operation.

### **Q. How do you access the memory stored in device?**

**A.** NONIN's nVISION data management software and 1000MC serial cable will allow you to connect the 2500 to PC for downloading data and view data.

### **Q. What happens to the stored memory after 72 hours have been displayed?**

**A.** The memory in the device functions as an "endless loop." When the memory fills up, the unit begins overwriting the oldest data with the newest.

## 2500A FAQ's

### **Q. What is included in the 2500A package assembly?**

- A.** • 2500A Pulse Oximeter • 2500A Operator's Manual on CD  
• 8000AA Pulse Oximeter Sensor • 4 AA-Size Alkaline Batteries

### **Q. How do you know you're getting an accurate reading?**

**A.** The tricolor LED display (m) provides a visual indication of the pulse signal quality, while blinking at the corresponding pulse rate. This display changes colors to alert you to changes in pulse quality that may affect the readings: *green* indicates a good pulse signal, *yellow* indicates a marginal pulse signal, and *red* indicates an inadequate pulse signal. Allow the 2500A to stabilize, observing about 10 seconds of continuous green-colored pulse quality before relying on the displayed values. It is common for the displayed values to fluctuate slightly over a period of several seconds. If the pulse quality display blinks yellow or red, try another finger

### **Q. What kind of battery is used in the 2500A?**

**A.** Four AA Alkaline or you can purchase the 2500C-UNIV Rechargeable NiMH Battery System.

### **Q. How many hours of battery life can I expect?**

**A.** Four AA alkaline batteries provide the device with approximately 80 hours of continuous operation, while the Rechargeable NiMH Battery Pack provides approximately 40 hours of continuous operation.

### **Q. How do you access the memory stored in device?**

**A.** NONIN's nVISION data management software and 1000MC serial cable will allow you to connect the 2500A to PC for downloading data and view data.

### **Q. Does this unit alarm?**

**A. Yes-** User's may adjust the alarm limits for upper and lower Spo2 and Pulse Rate alarms and alarm volume.