

Take the Fingertip Oximeter Challenge

1. Can your oximeter pass a very tough entrance exam?



Proof: The Onyx goes way beyond CE Marked and FDA Cleared requirements – including the new International Standard – ISO9919:2005. We take it to the next level passing stringent military aeromedical certifications – no other oximeter can make this claim. Ask for proof! www.nonin.com/onyx



2. Can you trust the accuracy of your oximeter to care for your most challenging patients?



Proof: Nonin welcomes scrutiny of our independent testing results and they are available for all to see! Check out why the Onyx is the most trusted fingertip oximeter at www.nonin.com/onyx.



3. Does your oximeter specifically list accuracy for patients with low perfusion?



Proof: Anyone can make low perfusion claims, but ONLY the Onyx has claims cleared by the FDA. Check out our accuracy specifications at www.nonin.com/onyx. Demand the same from the imitators.



4. Does your oximeter work on ANY patient – ANYTIME?



Proof: Nonin's PureSAT® oximetry technology is proven to be the most accurate on dark skin tones www.nonin.com/onyx. Don't forget the Onyx is a one-size fits all – we have the largest range of finger sizes from 0.3 to 1.0 inch thickness.



5. Can your oximeter survive the "50 drops" test?



Proof: Nonin's quality standards require every Onyx be able to withstand 50 drops from 1 meter onto a hard surface. Try it with an Onyx imitator and be prepared to pick up the pieces! View www.nonin.com/onyx to view the drop test video.



6. Will your oximeter still work when exposed to fluids?



Proof: The clinical environment can get messy. The Onyx exceeds IPX-2 water ingress testing. Be skeptical of the imitators – check out www.nonin.com/onyx to view the shower test video.



7. Do you expect to be using your oximeter in 5 years? 10 years?



Proof: Over 13 years of product innovations have resulted in less than 1% warranty return rate. People are still using the original Onyx from 1995 – today. How do imitators compare?



8. Does your oximeter company stand by its product?



Proof: The Onyx is NOT a throw-away product. If something happens, it can be repaired and like new in record time. Nonin is a company that will be here tomorrow.



9. Would you rather talk to a computer or a real person?



Proof: At Nonin, you'll speak with a real person right away. The voice you hear represents the pride in our products. We go above and beyond to have satisfied customers.

