



NiMH Battery Charger

NONIN's NiMH Battery Charger provides a convenient, economical solution to meet your battery needs.



Now you can easily charge NiMH batteries for NONIN's 2500 PalmSAT®, 8500 and 9840 Series of hand-held pulse oximeters. Charge a minimum of two, or up to a maximum of eight AA or AAA NiMH batteries in any combination of pairs. Fully charge depleted batteries in 8 hours, and automatically maintain the batteries in a fully charged condition without the risk of damage from overcharging. LED's indicate charge status: red when charging, off when batteries are fully charged.

Standard package includes charger and 12 high performance 1800 mAh NiMH batteries allowing you to power one or more NONIN pulse oximeters while charging additional batteries. The international charger model includes international power cord and adapters.

Product Highlights

- ✓ Cost Effective - Power one or more NONIN pulse oximeters while charging additional batteries
- ✓ Powerful - Fully charges depleted batteries in 8 hours and automatically maintains each battery in a fully charged condition
- ✓ Flexible - Can charge any combination of AA or AAA batteries at a time (in pairs)
- ✓ Worry Free - Charges AA or AAA NiMH batteries without the risk of damage from overcharging

Ordering Information

	Model #	Part #
Domestic NiMH Charger	NiMHBC-NA	4182-000
International NiMH Charger	NiMHBC-UNIV	4218-000
1800 mAh NiMH Batteries	NiMHB	4217-001



www.nonin.com

Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN 55441-5443, USA
Phone: +1.763.553.9968 Toll Free: 1.800.356.8874
Fax: +1.763.553.7807
info@nonin.com

✓ Where People Make the Difference