



Nonin RespSense™ - Sleep Lab Integration

Sens DAC Cable Matrix

| Manufacturer | Brand(s) | Voltage/Pin | Nonin P/N | Mfr. Contact Info |
|--------------------------------|--|--------------------------------------|--------------------------------------|---|
| Cadwell Laboratories | Easy II, Easy III | 0-1V/3.5mm phono | 7489-001 (2 pin) | U.S. toll free 800.245.3001 Tel. +1 509.735.6481 |
| Carefusion/Sensormedics | SomnoStar | 0-1V/1/8" phono | 7489-001 (2 pin) | U.S. toll free 800.231.2466 Tel. +1 714.283.2228 Germany +49 (0) 931 4972-0 |
| CleveMed | Crystal Monitor 20-S | 0-1V/Direct connection | 7491-000 | U.S. toll free 877.253.8363 Tel. +1 216.791.6720 |
| CompuMedics | E-Series | 0-1V/3.5mm phono | 7489-001 (2 pin) | U.S. toll free 877.717.3975 Tel. +1 704.749.3200 Germany +49 7731 79 76 9-0 |
| Embla | S4500, N7000 | 0-1V/RJ9 0-5V/RJ9 | 8123-001 (2 pin) 8805-001 (2 pin) | U.S. toll free 888.662.7631 International +31 20 3460130 |
| Grass Technologies | Comet | 0-1V/1/8" phono | 7489-001 (2 pin) | U.S toll free 877.472.7779 Tel. +1 401.828.4000 |
| Natus | Bio-logic, XLTEK | 0-1V/3.5mm phono | 7489-001 (2 pin) | U.S. toll free 800.303.0306 Tel. +1 905.829.5300 |
| Natus | Stellate | 0-1V/BNC | 7486-000 | U.S. toll free 800.303.0306 Tel. +1 905.829.5300 |
| NeuroVirtual | SleepVirtual BWII Sleep Virtual BWIII | 0-1V/3.5mm female 0-1V/3.5mm male | 8634-001 7489-001 | U.S. toll free 888-6-SLEEP-8 Tel. +1 786-693-8200 |
| Nihon Kohden | Neurofax | 0-1V/2.5mm mono 0-1V/2.5mm stereo | 7491-000 7487-000 | U.S. toll free 800.325.0283 Tel. +1 949.580.1555 |
| Respironics | Alice 5 | 0-1V/3.5mm phono | 7489-001 (2 pin) | U.S. toll free 800.345.6443 |

* Please contact Nonin Technical Service (tel. 1.800.356.8874) regarding PSG systems not listed above



Nonin LifeSense® - Sleep Lab Integration Sens DAC Cable Matrix

| Manufacturer | Brand(s) | Voltage/Pin | Nonin P/N | Mfr. Contact Info |
|--------------------------------|----------------------|--------------------------------------|--|---|
| Cadwell Laboratories | Easy II, Easy III | 0-1V/3.5mm phono | 7489-000 | U.S. toll free 800.245.3001 Tel. +1 509.735.6481 |
| Carefusion/Sensormedics | SomnoStar | 0-1V/1/8" phono | 7489-000 | U.S. toll free 800.231.2466 Tel. +1 714.283.2228 Germany +49 (0) 931 4972-0 |
| CleveMed | Crystal Monitor 20-S | 0-1V/Direct connection | 7491-000 | U.S. toll free 877.253.8363 Tel. +1 216.791.6720 |
| CompuMedics | E-Series | 0-1V/3.5mm phono | 7489-000 | U.S. toll free 877.717.3975 Tel. +1 704.749.3200 Germany +49 7731 79 76 9-0 |
| Embla | S4500, N7000 | 0-1V/RJ9 0-1V/RJ9 0-5V/RJ9 | 7488-000 (4 pin) 8123-001 (2 pin) 8805-001 (2 pin) | U.S. toll free 888.662.7631 International +31 20 3460130 |
| Grass Technologies | Comet | 0-1V/1/8" phono | 7489-000 | U.S toll free 877.472.7779 Tel. +1 401.828.4000 |
| Natus | Bio-logic, XLTEK | 0-1V/3.5mm phono | 7489-000 | U.S. toll free 800.303.0306 Tel. +1 905.829.5300 |
| Natus | Stellate | 0-1V/BNC | 7486-000 | U.S. toll free 800.303.0306 Tel. +1 905.829.5300 |
| NeuroVirtual | SleepVirtual | 0-1V/3.5mm female | 8634-001 | U.S. toll free 888-6-SLEEP-8 Tel. +1 786-693-8200 |
| Nihon Kohden | Neurofax EEG-9100 | 0-1V/2.5mm mono 0-1V/2.5mm stereo | 7491-000 7487-000 | U.S. toll free 800.325.0283 Tel. +1 949.580.1555 |
| Respironics | Alice 5 | 0-1V/3.5mm phono | 7489-000 | U.S. toll free 800.345.6443 |

* Please contact Nonin Technical Service (tel. 1.800.356.8874) regarding PSG systems not listed above