



# Nonin RespSense™ II - Sleep Lab Integration PSG DAC Cable Matrix

Manufacturer /Brand(s)	Voltage/Pin	Nonin P/N
<b>Cadwell Laboratories</b> (Easy II, Easy III)	0-1V/3.5mm Mono	10534-002 3.5 Mono, PSG DAC, RespSense II
<b>Carefusion/Sensormedics</b> (SomnoStar)		
<b>CompuMedics</b> (E-Series)		
<b>Natus</b> (Bio-logic, XLTEK)		
<b>Respironics</b> (Alice 5, Alice 6)		
<b>NeuroVirtual</b> (SleepVirtual)		
<b>Embla</b> (S4500, N7000)	0-1V/RJ9	10534-003 Embla PSG DAC
<b>Nihon Kohden</b> (Neurofax)	0-1V/2.5mm Mono	10534-004 2.5 Mono PSG DAC



# Nonin LifeSense® II - Sleep Lab Integration PSG DAC Cable Matrix

Manufacturer /Brand(s)	Voltage/Pin	Nonin P/N
<b>Cadwell Laboratories</b> (Easy II, Easy III)	0-1V/3.5mm Mono	10534-001 3.5 Mono, PSG DAC
<b>Carefusion/Sensormedics</b> (SomnoStar)		
<b>CompuMedics</b> (E-Series)		
<b>Natus</b> (Bio-logic, XLTEK)		
<b>Respironics</b> (Alice 5, Alice 6)		
<b>NeuroVirtual</b> (SleepVirtual)		
<b>Embla</b> (S4500, N7000)	0-1V/RJ9	10534-003 Embla PSG DAC
<b>Nihon Kohden</b> (Neurofax)	0-1V/2.5mm Mono	10534-004 2.5 Mono PSG DAC