

TRUST WHAT'S IN YOUR HANDS

Nonin CO-Pilot™

Wireless Handheld Multi-Parameter System (H500)



Pulse Oximetry



Pulse Rate



Carboxyhemoglobin



Methemoglobin



Cerebral or Tissue
Oximetry

Emergencies can happen anywhere.

And unreliable measurements leave everyone exposed.




You need to be able to depend on the tools you use so you can focus on the lives that depend on you. Introducing the Nonin CO-Pilot™ Multi-Parameter Wireless Handheld System to reliably measure five parameters including SpO2, COHb, and rSO2 in even the most challenging conditions. **So you can trust what's in your hands.**

A First-of-its-Kind Wireless, Multi-Parameter Monitoring System

Designed For Demanding EMS Environments

Engineered to withstand challenging conditions faced in the field, the Nonin CO-Pilot™ System is a flexible, five-parameter system capable of measuring SpO2, Pulse Rate, COHb, MetHb and rSO2 all on one wireless handheld device.

Responsive Readings

-  Accepts Nonin SpO2, rSO2 and multi-parameter sensors
-  Indicated for SpO2 accuracy from 80% to 100% with CO up to 15%
-  COHb accuracy with patients in normoxia and hypoxia

Reliability in the Field

-  Lightweight, handheld design is easy to carry and simple to use
-  Nonin multi-parameter and rSO2 sensors deliver dependable performance
-  Rugged design meets strict IEC 60601-1-12 EMS standards

Features Designed to Give Accuracy and Reliability

Proven durable sensors can monitor multiple parameters at once

Pleth waveform displays simultaneously with the parameters measured using the Multi-Parameter Sensor (8330AA)

Wireless display and signal processor are pre-paired; no additional steps needed

System can also monitor rSO2 during resuscitation or ECMO

Full battery charge provides approximately 10 hours of runtime



Nonin CO-Pilot™ Display and Carrying Case

Wireless Signal Processor

8330AA Fingertip Sensor

SenSmart® 8204CA Sensor

Nonin Medical, Inc.

13700 1st Avenue North
Plymouth, MN • 55441-5443 • U.S.A.

Tel: +1.763.553.9968 1.800.356.8874

Fax: +1.763.577.5521

Email: info@nonin.com

Nonin Medical B.V.

Doctor Paul Janssenweg 150
5026 RH Tilburg • Netherlands

Tel: +31 (0)13 – 45 87 130

Email: infointl@nonin.com



nonin.com



CE 0123

Rx Only