

ACCURACY FOR PEACE OF MIND

Nonin Model 7500FO (Fiber Optic) Pulse Oximeter



Easy to use



Lightweight



Durable



Versatile

Rest Easy. Act Confidently.

The Nonin Model 7500FO Pulse Oximeter empowers clinicians with accurate oximetry information ¹— helping them act quickly and confidently in the moments that matter.

The 7500FO is specifically designed for MR environments. Due to the 7500FO's fiberoptic sensor and cable design, radio frequency burns and image artifacts are eliminated.

Lightweight, compact, and portable, the 7500FO is ideal for monitoring infant, pediatric, and adult patients in fixed-site or mobile MR facilities.

A Monitor for What Matters

Pulse Strength Bar Graph

Assess the patient's status at a glance with a color-coded bar graph and a saturation level-related audible tone at each pulse.

Pulse Quality Indicator

Detects patient movement or signal quality changes that may require sensor site stabilization.

Sensor Fault Indicator

Indicates when a sensor has become disconnected, misaligned, failed, or is incompatible with the monitor.



Accuracy you can act on



PureSAT® Technology – Intelligent pulse-by-pulse filtering provides accurate oximetry measurements.²



Convenient – Lightweight with large, bright, LED display offers clear readability.



Powerful Internal Memory – 70 hours of memory available for extended use.



Flexible – Can be used in fixed-site or mobile MR facilities as a stand-alone device or interface with other equipment to monitor infant, pediatric, and adult patients.



Visit nonin.com for more information.

Authorized Distributor:

Nonin Medical, Inc.

13700 1st Avenue North
Plymouth, MN • 55441-5443 • U.S.A.

Tel: +1.763.553.9968 1.800.356.8874
Fax: +1.763.577.5521
Email: info@nonin.com

Nonin Medical B.V.

Doctor Paul Janssenweg 150
5026 RH Tilburg • Netherlands

Tel: +31 (0)13-45 87 130
Email: infointl@nonin.com