

# EXTENDING CARE BEYOND CLINICAL SETTINGS

## Model 3230 Bluetooth® Low Energy

OEM External Wireless Pulse Oximeter



Fast Readings



Easy Integration



Long Battery Life

The Nonin Model 3230 is an OEM external wireless finger pulse oximeter featuring Bluetooth® Low Energy (LE) technology. This high-performing OEM pulse oximeter provides SpO<sub>2</sub> and pulse rate spot-check monitoring. It features Bluetooth LE connectivity for simplified pairing and vital information exchange over a secure wireless connection.

# Accuracy in a Wide Range of Patients



Delivers accurate measurements for diverse skin pigmentation and even patients with low perfusion<sup>1</sup>



Proprietary CorrectCheck™ technology gives feedback when the finger is not placed properly in the device



SmartPoint™ sync feature indicated when the oximeter is syncing with the hub; user can remove oximeter after sync is complete



SmartPoint technology determines when a high quality measurement is ready to be transmitted



PureSAT® SpO<sub>2</sub> technology allows for pulse-by-pulse readings that can capture events such as rapid desaturations



Provides pulse amplitude index number in the data stream—no need to calculate your own number

## Your Partner Of Choice

More than 30 years of experience supporting OEM partners

More than 400 OEM partners globally

Provides up to 2,200 spot checks on two (2) AAA batteries

Dedicated integration support

High durability and quality backed by a 2-year warranty

Model 3230 Wireless Pulse Oximeter



Large patient-facing display

Automatic On/Off

Learn more at: [nonin.com/products/3230/](https://nonin.com/products/3230/) >>>

1. Nonin Medical, Inc. Clinical Data on File.

Nonin.com

United States  
info@nonin.com  
+1-763-553-9968

International  
info@nonin.com  
+31 (0)13-45 87 130

