

WristOx₂® Model 3150

Competency Checklist

Competency Objectives

To demonstrate proficiency in the clinical applications of monitoring the patient, set-up, and proper operation of the WristOx₂ Pulse Oximeter.

Participant will be required to demonstrate the following:

- Device set up for patient use
- Attachment of sensor to the patient
- Monitoring the patient
- Common troubleshooting techniques



Quick Setup	Completed	Initials
Unpack the WristOx ₂ Model 3150 and install the batteries		
Connect WristOx ₂ Model 3150 to the PC using nVision® for configuration		
Demonstrate activation button functionality		
Demonstrate how to select device activation options (Sensor activation, spot checking, and programmed)		
Select a sensor type		
Monitoring the Patient	Completed	Initials
Attach sensor		
Place WristOx ₂ Model 3150 on patient		
Review and describe each indicator on the display		
Troubleshooting	Completed	Initials
Demonstrate how to get back display when SpO ₂ and pulse rate are not showing		
Demonstrate what to do when the device will not pair		
Demonstrate what to do when dashes continually display on the LCD		
Optional Features	Completed	Initials
Show all screen display modes: <ul style="list-style-type: none"> • Full Display • Partial Display • Memory Volume (MVI) Display 		
Describe what each error code means (01, 02, 03, and 06)		
Explain and demonstrate how to clean the device		
Describe how to store the device after use		

Resources



WristOx₂®
Model 3150
with USB



WristOx₂®
Model 3150
with Bluetooth®
Low Energy

Employee Signature: _____ Date: _____

Evaluator Signature: _____ Date: _____

For additional information, visit [Nonin.com](https://www.nonin.com)

