

FOCUS ON WHAT MATTERS MOST

WristOx₂[®] Model 3150 with Bluetooth[®] Low Energy
Wrist-Worn Pulse Oximeter



Wireless Connectivity



Compact Design



Long Battery Life

Featuring Nonin PureSAT[®] SpO₂ technology, the WristOx₂ 3150 with Bluetooth[®] Low Energy (BLE) delivers fast, accurate oxygen saturation and pulse rate readings in challenging conditions, including low perfusion and motion¹. The WristOx₂ 3150 with BLE is compatible with a wide range of Nonin PureLight[®] sensors and provides continuous, reliable SpO₂ monitoring, and data recording². Applications include oxygen qualification and sleep diagnostic screenings conducted at home, in the hospital, and in sleep labs.

For Proven Performance, Count on Nonin.



Memory Volume Indicator (MVI) Mode provides a graphic display, verifying that a sufficient amount of data has been recorded



Data can be transferred via USB connection to nVISION® or wirelessly to the NoninConnect app



1,080 hours of non-volatile memory at four-second sampling rate



Dependable design and construction, backed by a 3-year warranty



Exclusive sensors pair easily with the WristOx₂ 3150 oximeter for optimal patient comfort



Long lasting, 44-hour battery life operating on two AAA batteries

Experience Nonin Quality and Greater Affordability

Small size wears like a watch with large display screen and enhanced sensor connector design

Transfer data securely and wirelessly to the NoninConnect app*, available on the Apple Store or download via USB cable

Single unit—includes sensor of choice and wristband starter kit—includes sensor of choice, three wristbands, carry case, USB download cable, and nVISION® software

WristOx₂ Model 3150 with Bluetooth® Low Energy



WristOx₂ Accessories



8000S-WO2 Reusable Sensor

*The NoninConnect app isn't cleared as a medical device

1,2. Nonin Medical, Inc. Data on File.

Nonin.com

United States
info@nonin.com
+1-763-553-9968

International
info@nonin.com
+31 (0)13-45 87 130

