

# STAY CONNECTED WITH YOUR HEALTH

## TruO2<sup>®</sup> Model 3245

Wireless Fingertip Pulse Oximeter with Bluetooth<sup>®</sup> Low Energy



Fast Readings



High Performance



Bluetooth Connectivity

The TruO2<sup>®</sup> Fingertip Pulse Oximeter is made with the same medical-grade pulse oximetry technology used in professional Nonin pulse oximeters for patients, clinicians, and the military. Bluetooth wireless connectivity allows caregivers to monitor patient information remotely, and make care decisions based on accurate and reliable patient measurements.<sup>1</sup>

**TruO2 is a Health and Wellness Device—No Prescription Required**

©2023 Nonin Medical, Inc. All trademarks are the property of Nonin Medical, Inc. unless otherwise noted.



# Accuracy. Connectivity. And Proven Performance.<sup>1</sup>



Provides quick, high-quality SpO<sub>2</sub> and pulse rate readings from small to large digits



Save SpO<sub>2</sub> history using a compatible Apple device and share spot check readings with the Apple Health App



Nonin PureSAT® technology captures SpO<sub>2</sub> and pulse rate measurements quickly and consistently



Pairs with NoninConnect App, allowing users to send SpO<sub>2</sub> history by email using a compatible Apple device



Nonin CorrectCheck™ technology provides feedback via digital display if the finger is not placed correctly in device



Nonin SmartPoint™ capture algorithm determines when a high-quality measurement is ready to be wirelessly transmitted

## Health App Compatible

TruO2® Model 3245  
Fingertip

Wireless connectivity allows viewing of SpO<sub>2</sub> and pulse rate from up to 10m (32ft) away on iOS® devices

Power turns on when applied and off when removed

Provides up to 2,200 spot checks on two AAA batteries<sup>1</sup>

Rugged, reliable design and construction backed by a 2-year warranty



Bright LED display

Automatic On/Off

1. Nonin Medical, Inc. Data on File.

**Nonin.com**

**United States**  
info@nonin.com  
+1-763-553-9968

**International**  
info@nonin.com  
+31 (0)13-45 87 130

