

RESPONSE® - SLEEP LAB INTEGRATION SENS DAC CABLE MATRIX

Manufacturer	Brand(s)	Voltage/Pin	Nonin P/N	Mfr. Contact Info
Cadwell Laboratories	Easy II, Easy III	0-1V/3.5mm phono	7489-001 (2 pin)	U.S. toll free 800.245.3001 Tel. +1 509.735.6481
Carefusion/Sensormedics	SomnoStar	0-1V/1/8" phono	7489-001 (2 pin)	U.S. toll free 800.231.2466 Tel. +1 714.283.2228 Germany +49 (0) 931 4972-0
CleveMed	Crystal Monitor 20-S	0-1V/Direct connection	7491-000	U.S. toll free 877.253.8363 Tel. +1 216.791.6720
CompuMedics	E-Series	0-1V/3.5mm phono	7489-001 (2 pin)	U.S. toll free 877.717.3975 Tel. +1 704.749.3200 Germany +49 7731 79 76 9-0
Embla	S4500, N7000	0-1V/RJ9 0-5V/RJ9	8123-001 (2 pin) 8805-001 (2 pin)	U.S. toll free 888.662.7631 International +31 20 3460130
Grass Technologies	Comet	0-1V/1/8" phono	7489-001 (2 pin)	U.S toll free 877.472.7779 Tel. +1 401.828.4000
Natus	Bio-logic, XLTEK	0-1V/3.5mm phono	7489-001 (2 pin)	U.S. toll free 800.303.0306 Tel. +1 905.829.5300
Natus	Stellate	0-1V/BNC	7486-000	U.S. toll free 800.303.0306 Tel. +1 905.829.5300
NeuroVirtual	SleepVirtual BWII Sleep Virtual BWIII	0-1V/3.5mm female 0-1V/3.5mm male	8634-001 7489-001	U.S. toll free 888-6-SLEEP-8 Tel. +1 786-693-8200
Nihon Kohden	Neurofax	0-1V/2.5mm mono 0-1V/2.5mm stereo	7491-000 7487-000	U.S. toll free 800.325.0283 Tel. +1 949.580.1555
Philips-Respironics	Alice 5, Alice 6	0-1V/3.5mm phono	7489-001 (2 pin)	U.S. toll free 800.345.6443

* Please contact Nonin Technical Service (tel. 1.800.356.8874) regarding PSG systems not listed above

