

RESPONSE™ II - SLEEP LAB INTEGRATION

SENS DAC CABLE MATRIX

Manufacturer /Brand(s)	Voltage/Pin	Nonin P/N
Cadwell Laboratories (Easy II, Easy III)	0-1V/3.5mm Mono	10534-002 3.5 Mono, PSG DAC
Carefusion/Sensormedics (SomnoStar)		
CompuMedics (E-Series)		
Natus (Bio-logic, XLTEK)		
Respironics (Alice 5, Alice 6)		
NeuroVirtual (SleepVirtual)		
Embla (S4500, N7000)	0-1V/RJ9	10534-003 Embla PSG DAC
Nihon Kohden (Neurofax)	0-1V/2.5mm Mono	10534-004 2.5 Mono PSG DAC

* Please contact Nonin Technical Service (tel. 1.800.356.8874) regarding PSG systems not listed above

