



FDA-cleared performance at your fingertips

TruO2[®] OTC Model 3250
Fingertip Pulse Oximeter with Bluetooth[®]



Accuracy across all skin tones



Available over the counter



Cleared by the FDA

The TruO2[®] OTC is a medical-grade pulse oximetry device that delivers accurate blood oxygen and pulse rate readings across all skin tones.¹ It's a durable, FDA-cleared device that provides pulse oximetry monitoring, available over the counter—no prescription required.

The inclusive fingertip oximeter that performs accurately across patients of all skin tones.¹



Accurate Oximetry. Over the Counter.



Inclusive pulse oximetry

Nonin's historic commitment to accurate pulse oximetry across all skin tones continues with the FDA-cleared device that is available without a prescription.¹



Easy tracking and sharing

Using a compatible iOS® device, save and share spot check readings with the Apple Health App, or send history by email through the NoninConnect™ App.



Fast and accurate readings

Equipped with Nonin PureSAT® technology, TruO2 OTC captures accurate oxygen and pulse rate measurements quickly.²



No guesswork

Nonin CorrectCheck™ technology notifies user if the finger is not placed correctly in the device.

Key Features

User friendly with simple operation

Bluetooth® wireless connectivity allows viewing of SpO₂ and pulse rate from up to 10m (32 ft) away on iOS® devices

Power turns on when applied and off when finger is removed

Bright and easy-to-read LCD display with clear digits and symbols

Provides up to 2,200 spot checks on two AAA batteries³

Rugged construction backed by a 2-year warranty

Wireless/Bluetooth®

Automatic On/Off

Bright LED/
Easy-to-Read Display

Long Battery Life

Durable Design



TruO2® OTC Model 3250 Fingertip

1, 2, 3. Nonin Medical, Inc. Data on File.

Nonin.com

United States
info@nonin.com
+1-763-553-9968

International
info@nonin.com
+31 (0)13-45 87 130

